


PspP I (Sau96 I isoschizomer)



5' ...G▼GNCC...3'
3' ...CCNG▲G...5'

PspPI is a restriction enzyme purified from *Psychrobacter immobilis* TA137.

Catalogue No 126-1, 1500 U
 126-2, 3x1500 U

Concentration 10-12u/μl and 40-60u/μl*

*Add an H to cat.# to order the high concentration

Reagents supplied: 10x M and 10x K buffer.

Unit substrate: Lambda DNA.

Unit calculation assay conditions: 50 mM NaCl, 10 mM Tris-HCl (pH 7.9 @ 25°C), 10 mM MgCl₂, 1 mM dithiothreitol, 100 μg/ml bovine serum albumin and DNA. Incubate at 25°C.

Absence of contaminants: 80 units of PspP I do not produce any unspecific cleavage products after 16 hrs incubation with 1 μg of Lambda DNA at 25°C. After 50-fold overdigestion with PspP I, greater than 95% of the DNA fragments can be ligated and recut with this enzyme.

Storage buffer: 50 mM NaCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 μg/ml bovine serum albumin and 50% glycerol. Store at -20°C.

Heat inactivation: 55°C for 15 minutes..

Methylation Sensitivity:

dam methylation: Not sensitive

dcm methylation: Blocked by overlapping

CpG methylation: Blocked by overlapping

Note: Incubation at 37°C results in 60% activity.

Reference: Rina, M., Caufrier, F., Mavromatis, K., Markaki, M., Kokkinidis, M. and Bouriotis, V. (1997) Gene, 197, 353-360.

Percent Activity in MINOTECH Buffers

L	M	H	SH	A	K
50-75	100	50	25-50	10	100


General reaction mixture:

10U PspPI	1μl
10x M or K buffer *	2μl
DNA substrate	<1μg
Sterile ultrapure water	Up to 20 μl
Incubate for 15 min at 25°C	

*In the case of M buffer we recommend the addition of BSA to a final concentration of 100 μg/ml.

Frequency of Cutting

λ	Ad-2	Φx174	pUC18	M13mp18	pBR322
74	164	2	6	4	15



Lambda DNA 1.0 % agarose