

Bgl I...



5' ...GCCNNNN▼NGGC...3'
3' ...CGGN▲NNNNCCG...5'

BglI is a restriction enzyme purified from *Bacillus globigii* lacking BglII.

Catalogue No 105-1, 2000 U
 105-2, 3x2000 U

Concentration 10-12u/μl and 40-
 60u/μl*

*Add an H to cat.# to order the high concentration

Reagents supplied: 10x U_{BglI} and 10x K buffer

Unit substrate: Lambda DNA.

Unit calculation assay conditions: 100 mM KCl, 20 mM Tris-HCl (pH 8.5 @ 25°C), 10 mM MgCl₂, 0.04% Triton X-100, 100 μg/ml BSA. Incubate at 37°C.

Absence of contaminants: 100 units of BglI do not produce any unspecific cleavage products after 16 hrs incubation with 1 μg of λ DNA at 37°C. After 50-fold overdigestion with BglI, greater than 95% of the DNA fragments can be ligated and recut with this enzyme.

Storage buffer: 300 mM NaCl, 20 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 500 μg/ml BSA and 50% glycerol. Store at -20°C.

Heat inactivation: 65°C for 20 minutes.

Methylation Sensitivity:

dam methylation: Not sensitive
dcm methylation: Not sensitive

CpG methylation: Blocked by some combinations of overlapping

Percent Activity in MINOTECH Buffers

L	M	H	SH	A	K
10-25	75-100	75-100	75-100	50	100

General reaction mixture:

10U BglI 1μl
10x U_{BglI} or K buffer * 2μl
DNA substrate <1μg
Sterile ultrapure water Up to 20 μl
Incubate for 15 min at 37°C

*In the case of U_{BglI} buffer we recommend the addition of BSA to a final concentration of 100 μg/ml.

Frequency of Cutting

λ	Ad-2	Φx174	pUC18	M13mp18	pBR322
29	20	0	2	1	3



Lambda DNA 0.7 % agarose