

MINOTECH RT

Catalogue No **801-1(10KU)**

Concentration **200U/μl**

Store at -20°C.

Reagents supplied: 5x MINOTECH RT assay buffer

Source: Purified from an *E. coli* strain carrying a plasmid with M-MuLV reverse transcriptase gene.

Description: High purity reverse transcriptase suitable for first strand cDNA synthesis.

Unit definition: One unit is defined as the amount of enzyme required to incorporate 1 nmol of dTTP into acid-insoluble material in a total reaction volume of 50 μl in 10 minutes at 37°C using poly(rA)•oligo(dT)18 as template.

Quality Control Assays:

- **Functional Assay:** MINOTECH RT is tested for performance in First strand cDNA synthesis followed by PCR with Taq polymerase. The resulting 1.600 bp PCR product is visualized as a single band on an ethidium bromide-stained agarose gel. *MINOTECH RT has been successfully used for synthesis of DNA fragments up to 8.8kb size.*

- **Absence of contaminants:** Tested extensively for the absence of nicking, endo- and exodeoxyribonucleases and RNases

Guaranteed stability: MINOTECH RT is guaranteed to maintain stability until expiration date

Recommended First-Strand cDNA

Synthesis mixture:

- 1 μl of oligo(dT) 20 (50 μM); or 200–500 ng of oligo(dT) 12-18 ; or 50–250 ng of random primers; or 2 pmol of gene-specific primer
- 10 pg–5 μg total RNA or 10 pg–500 ng mRNA
- 1 μl 10 mM dNTP Mix (10 mM each dATP, dGTP, dCTP and dTTP at neutral pH)
- Sterile ultrapure water up to 13 μl.

1. Heat mixture to 65°C for 5 minutes

2. Incubate on ice for at least 1 minute

3. Brief centrifugation and add:

- 4 μl 5X MINOTECH RT assay buffer
- 1 μl 0.1 M DTT
- 1 μl RNase Inhibitor (40 units/μl).
1 μl of MINOTECH RT (~200 units/μl)*

* 400 U of MINOTECH RT can be added to increase yield (for the generation of cDNA >5kb).

4. Mix gently. If using random primers, incubate tube at 25°C for 5 minutes.

5. Incubate at 37-42°C for 30–60 minutes.

6. Heat inactivation step at 70°C for 15 minutes.

Optional (recommended for PCR targets >1kb).

Remove RNA complementary to the cDNA, by adding 2 units of *E. coli* RNase H and incubate at 37°C for 20 minutes.

Recommended PCR mixture:

10x MINOTECH Taq pol. buf.	5 μ l
10mM dNTP mix	1 μ l
25 μ M forward primer	1 μ l
25 μ M reverse primer	1 μ l
cDNA(from first-strand reaction)	2 μ l
MINOTECH Taq DNA pol. (5 u/ μ l)	0.25-0.5 μ l
Sterile ultrapure water	Up to 50 μ l

Recommended PCR conditions:

Initial denaturation	94°C, 2min	
25-35 PCR Cycles	Denature	94°C, 45sec
	Anneal*	45-68°C, 30sec
	Extend	72°C, 1min/kb
Final extension	72°C, 10min	
Hold	4°C, indefinitely	

*Anneal temperature depends on primer T_m