


Sau3A I



5' ...▼GATC...3'
3' ...CTAG▲...5'

Sau3AI is a restriction enzyme purified from *Streptomyces* species.

Catalogue No 147-1, 200 U
 147-2, 1000 U
 147-3, 5000 U

Concentration 10-12u/μl and 40-60u/μl*

*Add an H to cat.# to order the high concentration

Reagents supplied: 10x M and 10x K buffer

Unit substrate: Lambda DNA.

Unit calculation assay conditions: 50 mM NaCl, 10 mM Tris-HCl (pH 7.9 @ 25°C), 10 mM MgCl₂, 1mM dithiothreitol, 100 μg/ml BSA. Incubate at 37°C.

Absence of contaminants: 50 units of *Sau3A I* do not produce any unspecific cleavage products after 16 hrs incubation with 1 μg of λ DNA at 37°C. After 50-fold overdigestion with *Sau3A I*, greater than 95% of the DNA fragments can be ligated and recut with this enzyme.

Storage buffer: 50 mM KCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 μg/ml BSA and 50% glycerol. Store at -20°C.

Heat inactivation: 65°C for 20 minutes.

Methylation Sensitivity:

dam methylation: Not sensitive

dcm methylation: Not sensitive

CpG methylation: Blocked by overlapping

Percent Activity in MINOTECH Buffers

L	M	H	SH	A	K
50	100	50	<10	50	100


General reaction mixture:

10U Sau3AI 1μl
10x M or K buffer * 2μl
DNA substrate <1μg
Sterile ultrapure water Up to 20 μl
Incubate for 15 min at 37°C

*In the case of M buffer we recommend the addition of BSA to a final concentration of 100 μg/ml.

Frequency of Cutting

λ	Ad-2	Φx174	pUC18	M13mp18	pBR322
116	87	0	15	7	22



Lambda DNA 1.4 % agarose