

Taq DNA Polymerase

Catalogue No 203-1, 500 u
 203-2, 1000 u
 203-3, 2500 u

Concentration 5 u/μl

Reagents supplied: 10x Taq DNA polymerase Buffer (1.5ml) and/or 10x Taq DNA polymerase w/o MgCl₂ (1.5ml) and 25mM MgCl₂ (1.5ml).

Source: Purified from an *E. coli* strain carrying a plasmid with Taq DNA polymerase gene from *Thermus aquaticus* YT-1.

Description: Taq DNA Polymerase is a thermostable enzyme that catalyzes 5'→3' synthesis of DNA. The enzyme has no detectable 3' → 5' proofreading exonuclease activity, but possesses low 5' →3' exonuclease activity.

Unit definition: One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nmoles of dNTPs into acid insoluble material in 30 minutes at 72°C.

Reaction conditions: 1x Taq polymerase buffer [50 mM KCl, 10 mM Tris-HCl pH 8.5 @ 25°C, 1.5 mM MgCl₂, 0.1% Triton X-100].

Quality Control Assays:

- Standard DNA Polymerase Assay Conditions (not PCR conditions):
The polymerase activity is assayed in 10 mM Tris-HCl (pH 8.5), 50 mM KCl, 1.5 mM MgCl₂, 200 μM each of dATP, dGTP, dCTP, dTTP (a mix of unlabeled and [³H] dTTP) and 12.5μg activated calf thymus DNA, in a final volume of 50 μl.

- Functional Assay: Taq DNA Polymerase is tested for performance in the polymerase chain reaction (PCR) using 1.5 units of enzyme to amplify a 1730-bp region of the PspP I methyltransferase gene from 5 ng of bacterial genomic DNA. The resulting PCR product is visualized as a single band on an ethidium bromide-stained agarose gel.

- Absence of contaminants: Tested extensively for the absence of endo- and exodeoxyribonucleases.

Guaranteed stability: Taq DNA polymerase is guaranteed to maintain stability for six months from the date of shipment when stored as directed.

Storage Buffer: 100 mM NaCl, 50 mM Tris-HCl (pH 8.0 @ 25°C), 1 mM DTT, 0.1 mM EDTA, 1% Triton X-100 and 50% glycerol. Store at -20°C.

Recommended PCR mixture:

10x Taq pol. buf.	5 μl
10mM dNTP mix	1 μl
25μM forward primer	1 μl
25μM reverse primer	1 μl
Template DNA	1-500 ng
Taq DNA pol. (5 u/μl)	0.25-0.5 μl
Sterile ultrapure water	Up to 50 μl

Recommended PCR conditions:

Initial denaturation	94°C, 2min
25-35 PCR Cycles	Denature 94°C, 45sec
	Anneal* 45-68°C, 30sec
	Extend 72°C, 1min/kb
Final extension	72°C, 10min
Hold	4°C, indefinitely

*Anneal temperature depends on primer T_m

Licenses/Patents/Disclaimers: *The PCR process is disclosed and claimed in European patents EP0201184B1 and EP0200362B1 or United States patent US4965188A, which have expired along with all family members/equivalents in all territories. The remaining technical features of this product do not appear to be covered by any in-force patent. However, some specific product applications may fall within the scope of protection of granted patents, which may be in force in the United States, European Union and certain countries. Because purchase of this product does not include a license to perform any application subject to patent protection, it is the sole responsibility of the buyer to ensure that use of the product does not infringe the patent rights of third parties. This product is for research use only.*