

Troubleshooting Guide for the Electrophoresis of DNA Markers

Problem	Probable Causes	Recommended Solutions
Faint or no DNA bands	<ul style="list-style-type: none"> • Quantity of DNA • DNA was degraded • DNA was electrophoresed off the gel • DNA was denatured 	<p>Increase the amount of DNA. A low concentration may be due to volume added per well width. Detection of DNA in polyacrylamide is less sensitive than in agarose.</p> <p>Avoid nuclease contamination of the DNA markers.</p> <p>Electrophorese the gel for less time, at a lower voltage, or in a higher percentage gel.</p> <p>Do not heat DNA markers (except Lambda-derived markers) prior to electrophoresis. Dilute markers in TE or in a buffer containing 20mM NaCl.</p>
Missing DNA bands	<ul style="list-style-type: none"> • Small DNA bands were electrophoresed off the gel • DNA bands of similar molecular size were not resolved 	<p>Electrophorese the gel for less time, at lower voltage, or in a higher percentage gel.</p> <p>Increase the electrophoresis time and check the proper percentage gel for resolution.</p>
Smearred DNA bands	<ul style="list-style-type: none"> • DNA was degraded • DNA was denatured • To much DNA was loaded on the gel • Inproper electrophoresis conditions were used • DNA contained too much salt • DNA was contaminated by protein 	<p>Avoid nuclease contamination of DNA markers.</p> <p>Do not heat DNA markers (except Lambda-derived markers) prior to electrophoresis. Dilute markers in TE or in a buffer containing 20mM NaCl.</p> <p>Decrease the amount of DNA in the gel.</p> <p>Do not allow voltage to exceed ~20 V/cm. Maintain a temperature <30°C during electrophoresis. Check that the electrophoresis buffer used has sufficient buffering capacity.</p> <p>Remove excess salt before electrophoresis by ethanol precipitation.</p> <p>Remove proteins before electrophoresis by phenol extraction.</p>
Anomalous DNA band migration	<ul style="list-style-type: none"> • DNA was denatured • Inproper electrophoresis conditions were used • Lambda DNA fragments, the cos site reannealed 	<p>Do not heat DNA markers (except Lambda-derived markers) prior to electrophoresis. Dilute markers in TE or in a buffer containing 20mM NaCl.</p> <p>Do not allow voltage to exceed ~20 V/cm. Maintain a temperature <30°C during electrophoresis. Check that the electrophoresis buffer used has sufficient buffering capacity.</p> <p>Heat DNA at 65°C for 5 min before electrophoresis.</p>